



STRUCTURED BREAST IMPLANT

SIZE	EMPTY + INNER + OUTER = TOTAL VOLUME					
	EMPTY IMPLANT VOLUME	INNER LUMEN SALINE	OUTER LUMEN SALINE	TOTAL IMPLANT VOLUME	DIAMETER	PROJECTION
			"HIGH" - "100%"	"HIGH" - "100%"	"HIGH" - "100%"	"HIGH" - "100%"
210cc	30cc	120cc	60cc - 85cc	210cc - 235cc	9.8cm - 9.7cm	4.0cm - 4.4cm
240cc	33	142	65 - 95	240 - 270	10.3 - 10.2	4.1 - 4.7
270cc	35	165	70 - 105	270 - 305	10.7 - 10.6	4.3 - 4.9
300cc	37	188	75 - 115	300 - 340	10.9 - 10.8	4.5 - 5.1
335cc	52	188	95 - 135	335 - 375	11.4 - 11.3	4.6 - 5.1
370cc	56	214	100 - 145	370 - 415	11.8 - 11.7	4.7 - 5.3
405cc	60	235	110 - 160	405 - 455	12.0 - 11.9	4.8 - 5.4
440cc	64	261	115 - 170	440 - 495	12.5 - 12.4	5.0 - 5.6
475cc	68	287	120 - 180	475 - 535	12.7 - 12.6	5.2 - 5.8
515cc	72	318	125 - 190	515 - 580	13.2 - 13.0	5.3 - 5.8
555cc	76	344	135 - 205	555 - 625	13.4 - 13.3	5.4 - 6.1
595cc	94	346	155 - 230	595 - 670	13.8 - 13.7	5.4 - 6.1
635cc	102	373	160 - 235	635 - 710	14.0 - 13.8	5.6 - 6.4
675cc	110	405	160 - 240	675 - 755	14.2 - 14.1	5.8 - 6.4

NOTE: Dimensions are on a curved surface. Conversion is necessary to compare to other manufacturers' dimensions on a flat surface, e.g., add 0.3cm to diameter and subtract 0.3cm from projection for a 210cc implant, add 0.5cm to diameter and subtract 0.5cm from projection for a 405cc implant, add 0.7cm to diameter and subtract 0.7cm from projection for a 675cc implant.



FEEL BEAUTIFUL AND SECURE

No risk of silent rupture

No MRI scans needed

Lower rupture rate

Lower capsular contracture rate


Natural feel and look

Safety of saline for peace of mind

For age 18 and older

Youthful shape and contour

THE STRUCTURED BREAST IMPLANT



Deflation/Rupture (MRI cohort) FDA Trial Data

PRIMARY AUGMENTATION		
IDEAL IMPLANT	7 yr	1.8%
Allergan Gel	7 yr	7.4%
Mentor Gel	6 yr	3.4%
Sientra Gel	6 yr	5.4%

Kaplan-Meier Risk Rates of First Occurrence

Capsule Contracture Baker III/IV FDA Trial Data

PRIMARY AUGMENTATION		
IDEAL IMPLANT	7 yr	6.3%
Allergan Gel	7 yr	16.2%
Mentor Gel	6 yr	9.8%
Sientra Gel	6 yr	10.0%

Kaplan-Meier Risk Rates of First Occurrence

Refer to the FDA-approved IDEAL IMPLANT Patient Information Booklet at idealimplant.com for information about the risks and benefits associated with the IDEAL IMPLANT® Structured Breast Implant. Additional information may be obtained at idealimplant.com.