

# IDEAL IMPLANT® STRUCTURED BREAST IMPLANT

Moderate Profile at Min. – Full Profile at Max.

EMPTY + INNER + OUTER = TOTAL VOLUME						
SIZE	EMPTY IMPLANT VOLUME	INNER LUMEN SALINE	OUTER LUMEN SALINE	TOTAL IMPLANT VOLUME	DIAMETER	PROJECTION
			MIN. – MAX.	MIN. – MAX.	MIN. – MAX.	MIN. – MAX.
210cc	30cc	120cc	60 - 85cc	210 - 235cc	10.1 - 10.0cm	3.5 - 4.3cm
240cc	33	142	65 - 95	240 - 270	10.5 - 10.4	3.6 - 4.5
270cc	35	165	70 - 105	270 - 305	11.0 - 10.8	3.8 - 4.7
300cc	37	188	75 - 115	300 - 340	11.4 - 11.2	3.9 - 4.9
335cc	52	188	95 - 135	335 - 375	11.9 - 11.7	4.0 - 5.1
370cc	56	214	100 - 145	370 - 415	12.2 - 12.0	4.1 - 5.2
405cc	60	235	110 - 160	405 - 455	12.5 - 12.4	4.2 - 5.4
440cc	64	261	115 - 170	440 - 495	12.9 - 12.8	4.3 - 5.6
475cc	68	287	120 - 180	475 - 535	13.3 - 13.1	4.4 - 5.7
515cc	72	318	125 - 190	515 - 580	13.6 - 13.4	4.5 - 5.8
555cc	76	344	135 - 205	555 - 625	13.9 - 13.8	4.6 - 6.0
595cc	94	346	155 - 230	595 - 670	14.3 - 14.2	4.7 - 6.1
635cc	102	373	160 - 235	635 - 710	14.6 - 14.5	4.8 - 6.2
675cc	110	405	160 - 240	675 - 755	14.9 - 14.8	4.9 - 6.3

Measured on a flat surface